

**SEM SOP**  
**Short Version**

Sample preparation	<p>Step1: Set the sample on the sample holder. There are several sample holders depending on sample size and shape.</p> <p>Step 2: Put the carbon tape on one end of the sample if sample is nonconductive.</p> <p>Note: Be careful about sample size (to high sample can damage sensitive parts of the SEM)</p>
User login and placing the sample holder at the SEM sample unit	<p>Step1: Write basic data to the user book (Name, date, time, account number, etc)</p> <p>Step2: Log in to the computer (Ctr+Alt+1 for the first screen and Ctrl+Alt+2 for the second screen)</p> <p>Step3: Hit icon on the desktop for SEM software window starting</p> <p>Step4: Put the gloves on</p> <p>Step5: Push VENT button and hold for a few second. When the light stop flashing unit is ready.</p> <p>Step6: Open sample unit (pull manually to open)</p> <p>Step7: Set the prepared sample holder on its position at the sample unit. Note: There is flat side of the sample holder-should be set forward during sample holder positioning. Set sample holder by one hand if possible</p> <p>Step8: Check whether the sample holder is set correctly</p> <p>Step9: Close the sample unit (manually- push forward to close). Hold the sample unit by one hand and push PUMP button and hold for a few seconds until you hear the sound of pumping.</p> <p>Note: during setting the sample and closing the sample unit the HT button at SEM window is not active.</p> <p>Note: Do not live any data at the desktop. Make sure that you have subfolder to storage your data.</p>
Sample investigation	<p>Step1: When HT button is active start your work</p> <p>Step2: Set the sample to the right position to have good view of the area you want to observe. Note: there are three ways to set the sample at the right position. Manually, with the screws at the SEM sample unit, with the buttons at the command table-joystick and by software (input X and Y coordinates). There is only one way to set Z coordinate manually by the screw at the SEM sample unit. Note: Be very careful with Z direction moving.</p> <p>Step3: Hit the SCAN2 button at the working window</p> <p>Step4: After setting Focus the sample (focus by minimum magnification)</p>

	<p>Note: Focus manually with the buttons at the command table (Focus + coarse = fast focusing). Focus can be done by using software buttons at the working window.</p> <p>Step5: Set the spot size around 60 and Acc Volt. at 20keV</p> <p>Step6: Set the brightness and contrast Note: There three ways for setting, by using manual command at the command table, by buttons at the working window menu and automatic settings with ACB button.</p> <p>Step7: Increase magnification to set the parameters for sample observation Note: Always use higher magnification for setting parameter than one you want to work on. (for example, use 1000X if you want to work with 100X)</p> <p>Step8: Focus the image</p> <p>Step9: Use OL-WOBBLER command from tools menu to set objective aperture properly. Note: correct manually with screws at the SEM sample unit</p> <p>Step10: Correct stigmation with buttons from command cable Note: Hit SCAN1 and use the same buttons as for contrast and brightness</p> <p>Step11: Focus the image</p> <p>Step12: decrease the magnification to your working magnification</p> <p>Step13: Focus the image</p> <p>Step14: Check contrast and brightness if necessary</p> <p>Note: If you want to see picture at the small window and compare some pictures press right click and chose snap shot option or import picture from the folder</p>
Working distance setting	<p>Step1: Check the Z coordinate and working distance at the working window Note: working distance is important for EDS analysis. Optimum working distance for EDS analysis is 12-14 mm Note: BE CAREFUL WITD SETTING WORKING DISTANCE</p> <p>Step2: Focus the image</p> <p>Step3: Decrease the working distance (decrease Z value for few mm)</p> <p>Step4: Focus the image</p> <p>Step5: Check working distance</p> <p>Repeat these steps until working distance reach 12-14 mm</p>

	Note: BE CAREFUL
Data saving	Step1: Hit SCAN3 and FREEZ Step2: Hit save data from the menu Step3: Open your Subfolder Step4: Give the name to the picture and save it
Working with BSEI	Note: use this option if you have multiphase sample  Step1: From the menu SIGNAL at working window chose BSEI  Step2: Chose appropriate option Compo or Topo or Shadow, depending on what you want to investigate  Note: When you finish return SIGNAL at SEI
EDS Analysis	Note: Sample must be flat and polished for this type of analysis  Step1: Go to screen 2 (Ctrl + Alt + 2) Hit the icon for starting EDS software  Step3: open the sample image at the screen 2  Step4: Chose the point where you want to analyze the sample  Step5: Follow the diagram for completing action for EDS analysis  Step6: Chose the type of data from the template menu you want to have at the end report  Note: Deadtime must to be between 40 to 60 % for correct results
Shutting of procedure	Step1: Set the minimum magnification  Step2: Decrease Z at the minimum position  Step3: Focus the image  Step4: Hit the HT button  Step5: Put the gloves on  Step6: Hit the VENT button  Step7: Open the Sample unit  Step8: Take off the sample holder  Step9: Close the sample unit and PUMP  Step10: Close the software

	Step11: Put the sample and sample holder at the right place
NOTE	If you have any comment about SEM or problem write it at the user book